

PREempt[®]

One-Step Surface Cleaner
and Disinfectant Wipes
formulated with



Wipes

WHY CHOOSE PREempt™?

- **CLEANER** Effective disinfectant as well as a cleaner resulting in added confidence that disinfection can occur.
- **FASTER** Rapid and realistic contact times on a broad spectrum of difficult to kill pathogens.
- **RESPONSIBLE** Designed to be easier on surfaces, staff and occupants. Users are not exposed to VOCs (volatile organic compounds) and PPE** is not required as per health and safety rating.
- **SUSTAINABLE** The active ingredient, Hydrogen Peroxide, breaks down into water and oxygen leaving no active residues which helps to reduce environmental impact.

DIRECTIONS FOR USE

Disinfection of non-porous hard surfaces such as the exterior of fume hoods, tabletops and other environmental surfaces.

1. Heavily soiled surfaces require cleaning prior to disinfection.
2. Apply to surface with disposable wipe.
3. *Ensure surface remains wet for 3 minutes.*

For use as a disinfectant on hard, non-porous environmental surfaces of footbaths:

1. **Drain:** Drain the water from the unit and remove any visible debris.
1. **Clean:** Thoroughly clean all surfaces with PREempt™ Wipes, then rinse with water and drain.
1. **Apply:** Pull towelettes from dispenser and wipe hard non-porous surfaces with the wipe.
1. **To Disinfect:** All surfaces must remain wet for 3 minutes. Use additional towelettes if needed to ensure the correct wet contact time.
1. **Rinse:** Rinse all surfaces with fresh water and drain.

For use in spas, esthetic schools, hair and nail salons, tanning salons, tattoo parlors, massage therapy centres and other critical environments that require cleaning and surface disinfection. This includes work stations, counter tops, chairs, spa and salon equipment and other hard non-porous environmental surfaces.

AVAILABLE PRODUCTS

PREempt™ One-Step Surface Cleaner and Disinfectant Wipes

DIN: 02436671

SHELF LIFE: 2 Years

CONCENTRATION OF ACTIVE:
0.5% Hydrogen Peroxide

- **Format Size:** 6"x7" – 160ct Wipes
Product Code: 11221
Packaging Specifications: 12x1 Case
- **Format Size:** 7.5"x10" – 10ct Wipes (Soft Pack)
Product Code: 11223
Packaging Specifications: 1 Case x 6 inner Cartons x 20 Soft Packs

EFFECTIVE AGAINST A BROAD-SPECTRUM OF BACTERIA, VIRUSES AND FUNGI IN 3 MINUTES!

Virucidal: 3 Min.

- Proven effectiveness against the Poliovirus Type 1, Sabin strain type 1 which allows for a Broad-Spectrum Virucidal claim against most enveloped and non-enveloped viruses
- HIV-1 Human Immunodeficiency Virus (HIV)
- Human Coronavirus 229E
- Feline Calicivirus, F9 Strain, as a surrogate for Norwalk and Norwalk-like viruses
- This product has demonstrated effectiveness against Poliovirus and is expected to inactivate all Influenza A viruses including 2009 (H1N1) pandemic Influenza A virus

Bactericidal: 3 Min.

- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- Methicillin-resistant *Staphylococcus Aureus* (MRSA)
- Vancomycin-resistant *Enterococcus faecalis* (VRE)

Fungicidal: 3 Min.

- *Trichophyton mentagrophytes*

Tuberculocidal: 3 Min.

- *Mycobacterium terrae*

Broad-Spectrum Sanitizing: 30 Sec.

- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- *Klebsiella pneumoniae*
- *Escherichia coli* O157:H7
- *Campylobacter jejuni*



PREempt™ One-Step Surface Cleaner and Disinfectant Wipes are non-woven, meltblown polypropylene wipe material, ensuring even surface coverage and increasing the ability and efficiency of AHP® for fine particle removal.

Please refer to reference sheet for a complete list of pathogenic organisms, additional features, and claims.

PREempt™ branded products are manufactured by Virox Technologies Inc.

Accelerated Hydrogen Peroxide® (AHP®) and PREempt™ are registered trademarks of Virox Technologies Inc.



For more information on this or our other products, visit:

> virox.com

2770 Coventry Road, Oakville, ON L6H 6R1 CANADA
1-800-387-7578 | www.virox.com



Engineering Revolutionary Disinfectants for the War Against Microbes