

PREempt®

Sterilant and High Level Disinfectant for Medical Devices and other Instruments formulated

with **ACCELERATED®**
HYDROGEN PEROXIDE
Patented



CS 20

WHY CHOOSE PREempt™?

- **CLEANER** The oxidation action of Hydrogen Peroxide provides added confidence that bio-burden potentially left behind during cleaning will not affix to surface.
- **FASTER** As a Ready to Use solution, there is no need to mix extra chemicals or add an activator. Rapid sterilization and high level disinfection in 20 minutes saves time and money.
- **RESPONSIBLE** Requires no special ventilation, is aldehyde free, and has no volatile organic compounds (VOCs). It is a more responsible solution for disinfection without the use of harmful chemicals.
- **SUSTAINABLE** The active ingredient, Hydrogen Peroxide, breaks down into water and oxygen, reducing environmental impact. May be reused for up to 14 days.

20-MINUTE STERILANT AND HIGH LEVEL DISINFECTANT 14 DAY REUSABLE SOLUTION

Sporicidal: 20 Min.

- *Bacillus subtilis*
- *Clostridium sporogenes*

Tuberculocidal: 20 Min.

- *Mycobacterium terrae*

Fungicidal: 20 Min.

- *Trichophyton mentagrophytes*

PREempt™ CS 20 is a sterilant and high level disinfectant for processing heat sensitive semi-critical and critical devices for which heat sterilization is not suitable.

AVAILABLE PRODUCTS

PREempt™ CS 20

LICENCE NO. 103055

SHELF LIFE: 12 Months

CONCENTRATION OF ACTIVE:
7% Hydrogen Peroxide

- 1. Format Size:** 1 L Bottle
Product Code: 11401
Packaging Specifications: 12x1 L
- 2. Format Size:** 4 L Bottle
Product Code: 11405
Packaging Specifications: 4x4 L

AHP7 Test Strips

- 1. Format Size:** 50 Test Strips per Bottle
Product Code: AHP7-S
Packaging Specifications: 1x1

DIRECTIONS FOR USE

Semi-Critical Devices:


1. Thoroughly clean the instruments and equipment mechanically or manually with an appropriate cleaning solution.
2. Rinse and dry objects before immersion in PREempt™ CS 20 solution.
3. Place pre-cleaned, rinsed and dried instruments and devices in solution of undiluted PREempt™ CS 20.
4. Allow devices to soak for 20 minutes. Rinse with sterile or potable tap water.

Critical Devices:

1. Thoroughly clean the instruments and equipment mechanically or manually with an appropriate cleaning solution.
2. Rinse and dry objects before immersion in PREempt™ CS 20 solution.
3. Place pre-cleaned, rinsed and dried instruments and devices in solution of undiluted PREempt™ CS20.
4. Allow devices to soak for 20 minutes. Remove and rinse with sterile water.

AHP7 Test Strips:

- AHP7 Indicator Strips were designed specifically for use with PREempt™ CS 20 to indicate the correct concentration of active.



PREempt™ CS 20 is powered by the patented Accelerated Hydrogen Peroxide® (AHP®) technology. AHP is a synergistic blend of commonly used, safe ingredients that when combined with hydrogen peroxide dramatically increases its germicidal efficacy and cleaning performance. This ensures that high level disinfection and sterilization can be achieved in realistic contact times while offering unparalleled user safety.

*Please refer to reference sheet for a complete list of pathogenic organisms, additional features, and claims.

PREempt™ brand products are manufactured by Virox Technologies Inc.

Accelerated Hydrogen Peroxide® and PREempt™ are registered trademarks of Virox Technologies Inc.



For more information on this or our other products, visit:

> virox.com

2770 Coventry Road, Oakville, ON CANADA L6H 6R1
1-800-387-7578 | www.virox.com

