

# PEROXigard®

## Wipes

One-Step Disinfectant Cleaner & Deodorizer

FOR USE IN LIFE SCIENCES

This product is a one-step disinfectant cleaner and deodorant designed for general cleaning, disinfecting, and deodorizing of hard, non-porous surfaces. It cleans quickly by removing urine, fecal matter, blood, food residue, dead skin, dirt, and other common soils found in lab animal research facilities, life science laboratories, small and large animal facilities, vivariums, animal housing, transportation vehicles, fish facilities, breeding establishments and operating rooms.

This product is designed for use on hard, non-porous surfaces such as: chrome, glass, glazed ceramic, glazed porcelain, glazed tile, plastic surfaces, stainless steel, and vinyl associated with high touch surfaces including countertops, biosafety cabinets, animal equipment, handheld devices, computer screens, cages, carts, and handling equipment.

## Claims

### BACTERICIDAL IN 1 MINUTE:

*Acinetobacter baumannii* (Multi-Drug Resistant (MDR))  
*Bordetella bronchiseptica*  
*Bordetella pertussis*  
*Burkholderia cepacia*  
*Chlamydia psittaci*  
*Enterobacter aerogenes*  
*Enterobacter cloacae* (New Delhi metallo-β-lactamase 1 (NDM 1))  
*Enterococcus faecalis* (Vancomycin Resistant (VRE))  
*Enterococcus faecium*  
*Escherichia coli*  
*Escherichia coli* O157:H7  
*Escherichia coli* (Extended Spectrum Beta-lactamase (ESBL))  
*Escherichia coli* (New Delhi metallo-β-lactamase 1 (NDM-1))  
*Klebsiella pneumoniae*  
*Klebsiella pneumoniae* Carbapenem Resistant (CRE)  
*Klebsiella pneumoniae* (Extended Spectrum Beta-lactamase (ESBL))  
*Klebsiella pneumoniae* (New Delhi metallo-β-lactamase 1 (NDM-1))  
*Legionella pneumoniae*  
*Listeria monocytogenes*

*Pasteurella multocida*  
*Pseudomonas aeruginosa*  
*Pseudomonas aeruginosa* (Metallo-β-lactamase (MBL) positive)  
*Salmonella enterica*  
*Serratia marcescens*  
*Shigella dysenteriae*  
*Staphylococcus aureus*  
*Staphylococcus aureus* (Community Acquired Methicillin Resistant Genotype USA 400 (CA-MRSA 400))  
*Staphylococcus aureus* (Methicillin Resistant (MRSA))  
*Staphylococcus aureus* (Vancomycin Intermediate Resistant (VISA))  
*Staphylococcus aureus* (Vancomycin Resistant (VRSA))  
*Staphylococcus epidermidis* (Methicillin Resistant (MRSE))  
*Staphylococcus pseudintermedius*  
*Stenotrophomonas (Pseudomonas) maltophilia*  
*Streptococcus pneumoniae* (Penicillin Resistant (PRSP))  
*Streptococcus pyogenes*

## Reference Sheet

KEEP OUT OF REACH  
OF CHILDREN

Manufactured by:

**Virox® Technologies Inc.**  
 2770 Coventry Road  
 Oakville, ON Canada L6H 6R1  
 1-800-387-7578 | virox.com



USER SHOULD READ S.D.S.  
BEFORE USING THE PRODUCT  
SDS available online at [www.virox.com](http://www.virox.com)

### \*VIRUCIDAL IN 1 MINUTE:

Adenovirus Type 8  
 Avian Influenza A Virus (H3N2)  
 Avian Influenza A Virus (H7N9)  
 Avian Reovirus  
 Bovine Viral Diarrhea Virus  
 Canine Distemper Virus  
 Canine Influenza A virus  
 Cytomegalovirus  
 Duck Hepatitis B Virus  
 Equine Herpes Virus Type 1  
 Equine Influenza A Virus (H3N8) (EIV)  
 Feline Calicivirus  
 Feline Herpes Virus (Feline Rhinotracheitis)  
 Feline Infectious Peritonitis Virus (FIPV)  
 Feline Leukemia Virus  
 Feline Panleukopenia Virus  
 Hepatitis B Virus  
 Hepatitis C Virus  
 Herpes Simplex Virus Type 1  
 Herpes Simplex Virus Type 2  
 HIV Type 1 (Strain HTLV III<sub>b</sub>)  
 Human Coronavirus  
 Human Coronavirus SARS-CoV-2 (the virus that may cause COVID-19<sup>2</sup>)  
 Influenza A Virus (Strain A/Hong Kong/8/68)  
 Influenza B Virus  
 Marek's Disease Virus  
 Minute Virus of Mice (Murine Parvovirus)  
 Mouse Hepatitis Virus  
 Murine Norovirus  
 Newcastle Disease Virus  
 Norovirus  
 Pandemic Influenza A Virus (H1N1)  
 Parainfluenza Virus Type 3  
 Poliovirus Type 1 (Strain Chat)  
 Porcine Epidemic Diarrhea Virus  
 Porcine Respiratory and Reproductive Syndrome Virus  
 Porcine Rotavirus  
 Pseudorabies Virus

Respiratory Syncytial Virus (RSV)  
 Rhinovirus Type-37  
 Rotavirus (WA strain)  
 Vaccinia Virus  
 Vesicular Stomatitis Virus

Meets surface disinfection recommendations under the OSHA Bloodborne Pathogen Standard for \*HBV, \*HCV & \*HIV-1

### TUBERCULOCIDAL IN 1 MINUTE:

*Mycobacterium bovis* (BCG)

### FUNGICIDAL IN 1 MINUTE:

*Candida albicans*  
*Microsporium canis*  
*Trichophyton interdigitale* (formerly known as *Trichophyton mentagrophytes*)

### NON-FOOD CONTACT SANITIZER IN 10 SECONDS:

*Acinetobacter baumannii*  
*Campylobacter jejuni*  
*Enterobacter aerogenes*  
*Enterococcus faecalis* (Vancomycin Resistant (VRE))  
*Escherichia coli* O157:H7  
*Listeria monocytogenes*  
*Pseudomonas aeruginosa*  
*Salmonella enterica*  
*Staphylococcus aureus*  
*Staphylococcus aureus* (Methicillin Resistant (MRSA))

### (MOLD/MILDEW) MILDEWCIDAL ACTIVITY IN 4 MINUTES:

*Aspergillus brasiliensis* (formerly known as *Aspergillus niger*)

### (MOLD/MILDEW) MILDEWSTATIC ACTIVITY IN 4 MINUTES:

*Aspergillus brasiliensis* (formerly known as *Aspergillus niger*)

# Directions For Use

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

**Not for personal cleansing. This is not a baby wipe.  
Do not flush in toilet.**

## For Use as a One-Step Cleaner/Disinfectant:

1. Pre-clean visibly soiled areas.
2. Pull towelette from canister and wipe hard, non-porous surfaces.
3. Allow surface to remain visibly wet for 1 minute.
4. Allow to air dry.

## For Use to Clean and Disinfect Life Science Laboratory Surfaces, Instruments, and Equipment:

1. Pre-clean visibly soiled areas.
2. Pull towelette from canister and wipe hard non-porous surfaces.
3. Allow surface to remain visibly wet for 1 minute.

## For Use in Facilities used for Temporary Confinement of Animals:

Use this product to clean and disinfect hard non-porous surfaces associated with facilities used for the temporary confinement of animals. Sites may include, but are not limited to, barns, sheds, stables, pens, cages, and associated access alleys or walkways. The product may also be used to clean and disinfect equipment related to the maintenance of animals found at fairs, exhibitions, animal auction yards, animal show/boarding facilities, or other similar agricultural facilities designed for temporary housing of animals.

1. Pre-clean visibly soiled surfaces and rinse with water.
2. Pull towelette from canister and wipe all hard, non-porous surfaces with this product.
3. Allow product to remain in contact with surfaces for 1 minute.

## For Use as a Cleaner/Disinfectant in Animal Housing Facilities:

For cleaning and disinfecting the following hard, non-porous surfaces: feeding and watering equipment, utensils, instruments, cages, kennels, etc.

1. Remove all animals and feeds from premises, animal transportation vehicles, crates, etc.
2. Remove all litter, droppings and fecal matter from floors, walls and surfaces of facilities occupied by animals.
3. Empty all feeding or watering appliances.
4. Thoroughly clean all surfaces with soap and detergent and rinse with potable water. Pull towelette from canister and wipe hard, non-porous surfaces.
5. All surfaces must remain visibly wet for 1 minute.
6. Allow to air dry.
7. Thoroughly scrub all treated feed racks, troughs, automatic feeders, fountains, and waterers with soap or detergent, and rinse with potable water before reuse.

## For Cleaning and Disinfection of Vehicle Cabins:

1. Pre-clean visibly soiled areas.
2. Pull towelette from canister and wipe hard, non-porous environmental surfaces.
3. Allow surface to remain visibly wet for 1 minute.
4. Allow to air dry.

## To Kill Mold and Mildew:

1. Pre-clean surfaces.
2. Pull towelette from canister and wipe product onto hard, non-porous surfaces until completely wet.
3. Allow the surface to remain visibly wet for 4 minutes.
4. Allow to air dry.

## To Control Mold and Mildew:

1. Pull towelette from canister and wipe product onto hard, non-porous surfaces until completely wet.
2. Allow to air dry.
3. Repeat application weekly or when growth reappears.

## To Disinfect Non-Critical<sup>¶</sup>, Pre-cleaned Instruments:

1. Instruments must be thoroughly pre-cleaned to remove excess organic debris, rinsed, and rough dried.
2. Clean and rinse lumens of hollow instruments before wiping with this product.
3. Wipe all hard non-porous surfaces of instruments with this product until thoroughly wet.
4. Let stand for 1 minute.
5. Allow to air dry.

## For use as a Non-Food Contact Sanitizer:

1. Pre-clean visibly soiled areas.
2. Pull towelette from canister and wipe hard, non-porous environmental surfaces.
3. Allow the surface to remain visibly wet for 10 seconds.
4. Allow to air dry.

## For Use as a Cleaner and/or Deodorizer:

1. Pull towelette from canister and wipe hard, non-porous environmental surfaces.
2. Allow to air dry.

**\*KILLS HIV-1, HCV AND HBV ON PRE-CLEANED ENVIRONMENTAL SURFACES / OBJECTS PREVIOUSLY SOILED WITH BLOOD / BODY FLUIDS** in health care settings and other settings in which there is an expected likelihood of soiling of hard non-porous surfaces / objects with blood / body fluids, and in which the surfaces / objects likely to be soiled with blood / body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1), Hepatitis C Virus (HCV) or Hepatitis B Virus (HBV).

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HCV OR HBV OF SURFACES / OBJECTS SOILED WITH BLOOD / BODY FLUIDS.**

**Personal Protection:** Disposable latex or vinyl gloves, gowns, face masks, and eye coverings must be worn during all cleaning of body fluids, blood, and decontamination procedures.

**Cleaning Procedures:** Blood / body fluids must be thoroughly cleaned from surfaces / objects before application of this product.

**Contact Time:** Allow surface to remain wet for 1 minute to kill HIV-1, HCV, HBV and all other organisms cited on label.

**Disposal of Infectious Material:** Cleaning materials used that may contain blood / body fluids should be autoclaved and disposed of according to Federal, State, and local regulations for infectious waste disposal.

\* on hard, non-porous surfaces.

¶ Non-Critical medical devices are items that come in contact only with intact skin.

\*Accelerated Hydrogen Peroxide® and design are trademarks of Diversey, Inc. | "Accelerated" does not imply increased product efficacy. | For relevant patents, visit: [virox.com/patents](http://virox.com/patents)

*Si no puede leer en inglés, pregunte a su supervisor sobre las instrucciones de uso apropiadas antes de trabajar con este producto. (If you cannot read English, ask your supervisor to explain the appropriate directions for use before working with this product.)*