



# PREvention®

Wipes

## One-Step<sup>†</sup> Disinfectant Cleaner FOR USE IN SPAS, SALONS & CLINICS

## Reference Sheet: Efficacy Claims & Use Directions



**KEEP OUT OF REACH  
OF CHILDREN**

READ S.D.S. BEFORE USING THE PRODUCT  
S.D.S. available online at [www.virox.com](http://www.virox.com)

Manufactured by:

**Virox® Technologies Inc.**  
2770 Coventry Road  
Oakville, ON CANADA L6H 6R1  
1-800-387-7578 | [virox.com](http://virox.com)

Meets surface disinfection recommendations under the OSHA Bloodborne Pathogen Standard for \*HBV, \*HCV & \*HIV-1  
Bactericidal, \*Virucidal, Fungicidal, and Tuberculocidal

### For Use In Spas, Salons & Clinics

Prevention™ Wipes are designed for cleaning and disinfection of hard non-porous surfaces in spas, nail salons, hair salons, health clubs, massage therapy clinics, medesthetic clinics, tanning facilities, tattoo parlors, and wellness centers. It cleans quickly by removing body oils, dead tissue, soil and all other buildups or organic matter found on hard non-porous surfaces.

It is designed for daily use on hard, non-porous environmental surfaces such as spa equipment, countertops, massage tables, workstations, facial tables, esthetic tools and equipment and non-circulating footbaths.

#### \*VIRUCIDAL: 1 MINUTE

Adenovirus Type 8  
Avian Influenza A Virus (H7N9)  
Avian Influenza A Virus (H3N2)  
Cytomegalovirus  
Feline Calicivirus  
Hepatitis B Virus (HBV)  
Hepatitis C Virus (HCV)  
Herpes Simplex Virus Type 1  
Herpes Simplex Virus Type 2  
HIV Type 1 (Strain HTLV III<sub>B</sub>)  
Human Coronavirus  
Human Coronavirus SARS-CoV-2  
(the virus that may cause COVID-19\*)  
Influenza A Virus (Strain A/Hong Kong/8/68)  
Influenza B Virus  
Norovirus (Feline Calicivirus as the surrogate)  
Pandemic Influenza A Virus (H1N1)  
Parainfluenza Virus Type 3  
Poliovirus Type 1 (Strain Chat)  
Respiratory Syncytial Virus (RSV)  
Rhinovirus Type-37  
Rotavirus (WA strain)  
Vaccinia Virus

#### BACTERICIDAL: 1 MINUTE

*Acinetobacter baumannii* (Multi-Drug Resistant (MDR))  
*Bordetella bronchiseptica*  
*Bordetella pertussis*  
*Burkholderia cepacia*

*Chlamydia psittaci* (Strain 6 BC)  
*Enterobacter aerogenes*  
*Enterobacter cloacae* (New Delhi metallo-β-lactamase 1 positive (NDM -1))  
*Enterococcus faecalis* (Vancomycin Resistant (VRE))  
*Enterococcus faecium*  
*Escherichia coli*  
*Escherichia coli* O157:H7  
*Escherichia coli* (Extended Spectrum Beta-lactamase (ESBL) positive)  
*Escherichia coli* (New Delhi metallo-β-lactamase 1 (NDM-1) positive)  
*Klebsiella pneumoniae*  
*Klebsiella pneumoniae* (Carbapenem Resistant (CRE))  
*Klebsiella pneumoniae* (Extended Spectrum Beta-lactamase (ESBL) positive)  
*Klebsiella pneumoniae* (New Delhi metallo-β-lactamase 1 (NDM-1) positive)  
*Legionella pneumoniae*  
*Listeria monocytogenes*  
*Pasteurella multocida*  
*Pseudomonas aeruginosa*  
*Pseudomonas aeruginosa* – (Metallo-beta lactamase (MBL) positive)  
*Salmonella enterica* (formerly known as *Salmonella choleraesuis*)  
*Serratia marcescens*  
*Shigella dysenteriae*  
*Staphylococcus aureus*  
*Staphylococcus aureus* (Community Acquired Methicillin Resistant Genotype USA 400)  
*Staphylococcus aureus* (Methicillin Resistant (MRSA))  
*Staphylococcus aureus* (Vancomycin Intermediate Resistant (VISA))  
*Staphylococcus aureus* (Vancomycin Resistant (VRSA)) (VRS1)  
*Staphylococcus epidermidis* (Methicillin Resistant (MRSE))  
*Staphylococcus pseudintermedius*  
*Stenotrophomonas (Pseudomonas) maltophilia*  
*Streptococcus pneumoniae* (Penicillin Resistant (PRSP))  
*Streptococcus pyogenes*

#### TUBERCULOCIDAL: 1 MINUTE

*Mycobacterium bovis* (TB) (BCG)

#### FUNGICIDAL: 1 MINUTE

*Candida albicans*  
*Microsporium canis*  
*Trichophyton interdigitale* (formerly known as *Trichophyton mentagrophytes*) (athlete's foot fungus)

## NON-FOOD CONTACT SANITIZER: 10 SECONDS

*Acinetobacter baumannii*

*Campylobacter jejuni*

*Enterobacter aerogenes*

*Enterococcus faecalis* (Vancomycin Resistant (VRE))

*Escherichia coli* O157:H7

*Listeria monocytogenes*

*Pseudomonas aeruginosa*

*Salmonella enterica* (formerly known as

*Salmonella choleraesuis*)

*Staphylococcus aureus*

*Staphylococcus aureus* (Methicillin Resistant)(MRSA)

## (MOLD/MILDEW) MILDEWCIDAL ACTIVITY: 4 MINUTES

*Aspergillus brasiliensis* (formerly known as *Aspergillus niger*)

## (MOLD/MILDEW) MILDEWSTATIC ACTIVITY: 4 MINUTES

*Aspergillus brasiliensis* (formerly known as *Aspergillus niger*)

## DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

Not for personal cleansing. This is not a baby wipe.  
Do not flush in toilet.

### For Use as a One-Step Cleaner/Disinfectant:

1. Pre-clean visibly soiled areas.
2. Remove towelette from canister and wipe hard, non-porous surfaces.
3. Allow surfaces to remain visibly wet for 1 minute.
4. Allow to air dry.

### For Use as a Non-Food Contact Sanitizer:

1. Pre-clean visibly soiled areas.
2. Remove towelette from canister and wipe hard, non-porous environmental surfaces.
3. Allow the surface to remain visibly wet for 10 seconds.
4. Allow to air dry.

### For Use as a Cleaner and/or Deodorizer:

1. Remove towelette from canister and wipe hard, non-porous environmental surfaces.
2. Allow to air dry.

### For Use as a Disinfectant, on Hard, Non-Porous Surfaces of Non-Circulating Footbaths:

1. **Drain:** Drain the water from the basin and remove any visible debris.
2. **Clean:** Thoroughly clean all hard, non-porous surfaces with this product, then rinse with water and drain.
3. **Apply:** Wipe all hard non-porous surfaces. Allow surfaces to remain visibly wet for 1 minute.
4. **Rinse** Surfaces with fresh water.

### To Disinfect Non-Critical<sup>¶</sup>, Pre-cleaned Instruments:

1. Instruments must be thoroughly pre-cleaned to remove excess organic debris, rinsed, and rough dried.
2. Wipe all hard non-porous surfaces of instruments with this product until thoroughly wet.
3. Let stand for 1 minute.
4. Allow to air dry.

### To Kill Mold and Mildew:

1. Pre-clean surfaces.
2. Remove towelette from canister and wipe product onto hard, non-porous surfaces until completely wet.
3. Allow the surfaces to remain visibly wet for 4 minutes.
4. Allow to air dry.

### To Control Mold and Mildew:

1. Remove towelette from canister and wipe product onto hard, non-porous surfaces until completely wet.
2. Allow to air dry.
3. Repeat application weekly or when growth reappears.

**\*KILLS HIV-1, HCV AND HBV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings and other settings in which there is an expected likelihood of soiling of hard, non-porous surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis C Virus (HCV) or Hepatitis B Virus (HBV).

### SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV OR HCV ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:

**Personal Protection:** Disposable latex or vinyl gloves, gowns, face masks, and eye coverings must be worn during all cleaning of body fluids, blood, and decontamination procedures.

**Cleaning Procedures:** Blood/body fluids must be thoroughly cleaned from surfaces and objects before application of this product.

**Contact Time:** Allow surface to remain visibly wet for 1 minute to kill HIV-1, HCV and HBV.

**Disposal of Infectious Material:** Cleaning materials used that may contain blood/body fluids must be autoclaved and disposed of according to Federal, State, and local regulations for infectious waste disposal.

## STORAGE AND DISPOSAL

**PESTICIDE STORAGE:** Store this product in a cool dry area, away from direct sunlight and heat. Do not contaminate water, food or feed by storage and disposal.

**CONTAINER HANDLING:** Dispose of wipe in trash after use. Do not flush. Tightly close lid between uses to retain moisture. Do not reuse or refill this container/dispenser. Offer empty container/dispenser for recycling. If recycling is not available, discard container/dispenser in trash.

<sup>¶</sup> Non-Critical medical devices are items that come in contact only with intact skin.

<sup>†</sup> When used according to disinfection use directions.

<sup>‡</sup> on hard, non-porous surfaces

<sup>©</sup> Accelerated Hydrogen Peroxide<sup>®</sup> and design are trademarks of Diversey, Inc. | For relevant patents, visit: [virox.com/patents](http://virox.com/patents)